140 mm

36 mm

18 mm

400 mm

12 mm

VisuXL[®] gel

Ophthalmic gel

EN Instructions for use

VisuXL[®] Gel is a preservative free ophthalmic gel with lubricating and antioxidant components: Coenzyme Q10, Vitamin E TPGS and cross-linked sodium carboxymethylcellulose.

VisuXL[®] Gel is indicated for both day and night use. Cross-linked sodium carboxymethylcellulose is a biosynthetic substance, safe and inert, effective and long lasting.

Its cross-linking and the lack of specific enzymes able to degrade the molecule at the level of the ocular surface, stabilises the polymer giving it a longer half-life compared with the traditional hyaluronic acid.

Furthermore, the cross-linked sodium carboxymethylcellulose is a hydrogel which in its native state has specific characteristics including muco-mimetic properties, viscosupplementation, hydration and substance delivery.


The combination of cross-linked sodium carboxymethylcellulose and antioxidant/anti free-radical substances makes VisuXL[®] Gel a highly effective treatment in case of:

- Moderate to severe dry eye;
- Foreign body trauma;
- Alterations in the continuity of the corneal and conjunctival surfaces after refractive corneal surgery and corneal transplant. In order to alleviate the typical post-surgical symptoms, VisuXL[®] Gel acts as a cushion between the eyelid and the corneal and conjunctival epithelia;
- Alterations of ocular surface related to metabolic disorders (diabetes).

HOW TO USE

Instil 1 drop of VisuXL[®] Gel in the conjunctival sac of each eye 2 times a day or as directed by your ophthalmologist.

How to apply VisuXL[®] Gel

WARNINGS AND PRECAUTIONS

- Do not use this product if you are allergic or hypersensitive to any of the components.
- If you have any adverse reactions, discontinue its use and contact your ophthalmologist.
- The product can also be instilled while wearing contact lenses: it does not alter their characteristics, and can even hydrate soft contact lenses.
- Instil the product at least 5 minutes after the use of other drugs.
- Do not use this product after the expiry date printed on the box.
- The product remains sterile up to 180 days after first opening of the bottle.
- The expiration date refers to the product in an unopened, correctly stored package.
- Wash your hands before applying the product.
- Do not let the tip of the bottle touch your eyes.
- Do not use the bottle if damaged.
- Keep the product away from light and heat and out of the reach of children.
- Store the product at room temperature (15-25 °C). Ophthalmic use, do not swallow.

COMPOSITION

- 100 ml of VisuXL[®] Gel contains:
- Coenzyme Q10 100 mg;
 - Vitamin E TPGS (D-α-Tocopherol polyethylene glycol 1000 succinate) 500 mg;
 - Cross-linked sodium Carboxymethylcellulose 400 mg;
 - Poloxamer 407 9000 mg;
 - Disodium EDTA 100 mg;
 - Buffered isotonic solution q.s. to 100 ml.

PACKAGE

Bottle of 10 ml, preservative free in 3K System
Bottle of 5 ml, preservative free in 3K System
Bottle of 3 ml (free sample), preservative free in 3K System

DISTRIBUTOR

VISUFARMA B.V.
Amstelplein 1, 1096 HA, Amsterdam, Netherlands
info@visufarma.com

MANUFACTURER

VISUFARMA S.p.A.
Via A. Cadlolo 21, 00136 Rome, Italy
info@visufarma.it

DO NOT DISPOSE OF THE PRODUCT IN THE ENVIRONMENT AFTER USE

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DISTRIBUTEUR VISUFARMA B.V.
FABRICANT VISUFARMA S.p.A.